







Primo Passi Pump Data- Copyright Genuine Lactation LLC 2023

	Double Pumping Original Flanges	BB Version 1
Stim 1		
	50 Short- 280 Long	
	82 cpm	




65 Short- 280 Long

74 cpm




80 Short- 280 Long

70 cpm




Stim 4




105 Short- 280 Long

68 cpm




Stim 5



125 Short- 280 Long


51 cpm





Exp 1

70

76 cpm





Exp 2



120


57 cpm

Exp 3




160

45 cpm




Exp 4




190

33 cpm






Exp 5



220


28 cpm





Exp 6

250

22 cpm





Exp 7



285



18 cpm

Exp 8



325

15 cpm



Exp 9

350

14 cpm